EYE CARE PROFESSIONAL FITTING GUIDE


Jena Indoor lenses are the clearest and most precise lenses available for indoor use.

By using your specific facial measurement, frame parameters and prescription, our German custom design software will engineer personalized lenses with almost no distortion – perfect comfort.

INDOOR DESIGN SELECTIONS

PATIENT SELECTION

- Patients who need some intermediate vision.
- Patients who dislike or are new to progressive lenses.
- Patients who accept that they may need two pairs of eyeglasses*
*(*or will have to remove the Office glasses if they can see Far without a prescription. Or will have to have to wear contact lenses.)*

HANDLING FAR VISION REQUIREMENTS

For patient satisfaction, a discussion must be made to find out how they will handle their far vision needs

- Removing Office Lens (plano Far Rx)
- Contact lenses (office lens must be converted to accommodate contact lens over refraction).
- Separate set of glasses.
- Are they willing to accept some diminished distance? (non drivers only)

FITTING & DISPENSING

All Jena products utilize true backside freeform optics. The lens has been engineered so that some parameters are more critical than others.

Place frame on patient and adjust frame before measurements are taken.

- Mono PD (exact)
- Height to pupil (round down if unsure)
- Pantoscopic Tilt (approximate)
- Vertex Distance (approximate)
- Wrap angle (approximate)

Minimum fitting height should be no less than 17mm. The prescription on the lens will be manufactured to 0.01D.

DEVICE

Up to 1 meter (3'); Enhanced Readers.

Intermediate Rx = Distance Rx + (55% of ADD)

WORKSPACE

Up to 2 meters (6'); Popular with heavy computer users.

Intermediate Rx = Distance Rx + (38% of ADD)

STUDIO

Up to 4 meters (12'); For complete indoor vision.

Intermediate Rx = Distance Rx + (24% of ADD)

AVAILABILITY

● = Available with Zero Blue technology!

MATERIAL	POWER RANGE	ZERO BLUE
THINNER	1.50 - 5.00 to + 4.00 (+0.75 to +3.50)	
	Trivex - 6.00 to + 5.00 (+0.75 to +3.50)	
	1.56 - 6.00 to + 5.00 (+0.75 to +3.50)	●
	Poly - 10.00 to + 7.00 (+0.75 to +3.50)	●
	1.6 - 10.00 to + 7.00 (+0.75 to +3.50)	●
	1.67 - 13.00 to + 7.50 (+0.75 to +3.50)	●
	1.74 - 15.00 to + 8.00 (+0.75 to +3.50)	●

High Energy Blue blocking coatings are recommended or a high quality AR coating. Protect your patient's eyes and health with an upgrade to **Zero Blue** technology, the most effective high energy blue light protection available!

Included in the **Zero Blue** coating are an invisible high energy **Clear Blue™** reflective coating and an absorbing **Blue Barrier™** material.

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